Technovit® 2000 LC

DE Gebrauchsinformation Technovit 2000 LC System

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GB User instructions Technovit 2000 LC system

Application

Technovit 2000 LC is a light-curing resin for embedding of translucent ground sections within material testing. Technovit 2000 LC had been specifically developed for preparation of sensitive materials and micro components. The main application areas are: Semiconductor technique, microelectronics, medical equipment technique, optoelectronics and microsystem techniques. Technovit 2000 LC varnish is a light-curing varnish to avoid the development of an inhibition layer on ground sections embedded in Technovit 2000 LC. Application in the laboratory by trained specialists.

Material

Light-curing resin in the form of liquid based of polyfunctional methacrylates. Technovit 2000 LC and Technovit 2000 LC varnish are 1-component Materials.

Properties

The Technovit 2000 LC liquid is light-curing and is polymerized within the none-dangerous spectrum of the visible blue-light generated by the polymerization unit Technotray POWER of Kulzer. The arbitrarily long pot life allows for optimal positioning and infiltration of the sample (or removal of any air bubbles). If properly handleed gap-free embeddings will be the result. The cured Technovit 2000 LC can be optimally processed mechanically.

Advantages

- dry and non-sticky surface of the ground section embedded in Technovit 2000 LC
- no loss of material (one-component resin)
- light-curing thus almost unlimited handling time
- optimum translucency
- exceptional low gap formation
- no air entrapplings
- alcohol and acid resistant
- usable under vacuum, no gas evaporation (can be used under REM)
- odor-free
- low polymerization temperature

Handling

Surfaces exposed to Technovit 2000 LC must be clean and degreased for embedding. For this degreasing solvents such as light benzine or ester are recommended. The Kulzer embedding molds are used for embedding; silicone molds are NOT suitable.

The samples can be positioned and fixed with the light-curing Technovit 2000 LC fixing paste.

Curing

Polymerization is possible in block or in layers. Larger specimens should be embedded in layers. Polymerization takes place in the curing unit Technotray POWER. This unit is specifically tuned to Technovit 2000 LC.

Set the timer to 10 minutes and close the tray. After completion of the polymerisation cycle, pour on the Technovit 2000 LC thinly (max. 2 mm) so that the surface of the sample is completely wetted. The covering varnish is used for prevention of the dispersion layer on the top of the sample. Afterwards, start another polymerisation cycle of 10 minutes. After completion of the polymerisation time, remove the sample, let it cool and then divest it.

Working with Technovit 2000 Inside Cure

In order to harden Technovit 2000 LC in porous samples, shady places or lighting-inaccessible places (internal threads, pipe sections, capacitors, etc.), Technovit 2000 LC is mixed with the Inside Cure additive in a mixture of 1000 ml to 40 ml. For this purpose pour the contents of the original packing of Technovit Inside Cure (40 ml) into an original packing of Technovit 2000 LC (1000 ml). Tightly seal the Technovit 2000 LC bottle and shake well.

As this mixture has a shorter shelf life than Technovit 2000 LC please place the "Inside Cure mixture" label accompanying the packing of the Technovit Inside Cure on to the Technovit 2000 LC bottle and note the life from the beginning of the mixture as the shelf life + 1 year.

Technovit 2000 LC can now be used as usual but it is not recommended to work in temperature-time programmes as then the polymerisation cannot take place in a regulated manner. In this case, it must not be processed in layers.

Delivery units

Technovit 2000 LC	– 1000 ml Liquid
Technovit 2000 LC varnish	– 100 ml Liquid
Technovit 2000 LC fixing paste	– 1 x 4 g
Technovit Inside Cure	– 1 x 40 ml Liquid
Technotray POWER light unit	

Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.

www.kulzer-technik.com

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery. If the assumption is justified that there are quality defects caused by the manufacturer, please notify us immediately.

* For more information please visit www.kulzer-technik.com

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© = marchio registrato, Kulzer GmbH, Hanau