- (DE) Gebrauchsinformation
- (GB) User instructions
- (FR) Mode d'emploi
- (IT) Informazioni per l'uso
- (ES) Modo de empleo

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## (GB) User instructions Technovit 2-Component Resins

The fast and coldcuring Technovit resins in form of powder and liquid are state of the art in their field of application.

### Advice for materialographic testings:

The specimen must be clean and free of grease before embedding. For cleaning ethanole or ben zine is recommended. Depending on the specimen pour a thin layer of Technovit resin into the embedding mould first. Then put in the specimen and fill up the mould until the margin. Easy shaped specimen can be placed onto the bottom of the embedding mould and doused with Technovit resin directly. Let flow out the resin in a thin stream to allow airbubbles to pass off. Polymerized resins can be treated with all usual mechanical methods. Especial for materialographic applications the grinding and polishing products from the Kulzer series achieve best results. Application in the laboratory by trained specialists.

### Mixing

**A homogeneous mixture of the resin is the precondition to achieve optimal properties!** Place liquid, then put in the half of the powder and mix thoroughly. Add the 2nd half of powder after and after and mix for 40 seconds. Let the dough set down for 10–15 seconds and then stir again for at least 20 seconds. Basic mixing ratio for 2-component-resins see table; ratio can be varied according to the requirements. Stir thoroughly with a spatula until the powder is completely soaked with the liquid. Avoid beating movement to include airbubbles in the resin dough. A few hard taps on the worktable will help to escape air out of the resin. For mixing Technovit hardpaper cups and wooden spatulas are recommended. The heat caused by polymerization can my lead to formation of pores. To avoid excess heat do not mix more than 30g of dough (20 g powder/10 g liquid). Carry off heat by cooling (in a vessel with cold water — do not immerse) is recommended. Larger specimen has to be prepared in several layers. Please be sure, that each layer is complete polymerized before applying the next one. For absolutely bubble free embeddings, use the Technomat pressure pot.

**Storage advice** Keep vessels always closed. Do not store above room temperature. After expiry date the materials should not be used anymore. Keep away from reach of children.

### Danger/Safety advices Please follow the instructions on the product packaging and safety data sheets.

[www.kulzer-technik.com](http://www.kulzer-technik.com) Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery. If the assumption is justified that there are quality defects caused by the manufacturer, please notify us immediately. Shelf life see product packaging. Please adhere to the temperatures specified in the processing instructions to ensure error-free use.

\* For more information please visit [kulzer-technik.com](http://kulzer-technik.com)

® = registered trademark, Kulzer GmbH, Hanau

66094619\_01 - 11/2022  
Technovit Kombi incl. Technovit 4021 - 5-sprachig

**Manufacturer:**  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau, Germany  
Made in Romania

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® = marchio registrato, Kulzer GmbH, Hanau
® = marca registrada., Kulzer GmbH, Hanau